LABORATORY PRODUCTS



ACROPARS Re

Cold Cure Rezin Acrylic

Acropars Re is designed for repairing or extending the full or partial prosthesis and is based on methyl methacrylate copolymer. It is comply with ISO 1567 and ADA No. 12.



- 1- Repairing denture prosthesis.
- 2- Extending denture prosthesis.
- 3- Relining of denture prosthesis.

Re Actionals Cold-core acrylic for respect The pass and the pass and

ADVANTAGE:

- 1- Enough time to repair or reline.
- 2- Very good adhetion to old prosthesis.
- 3- Good polish ability.
- 4- Good color stability.

PACKAGING:

Powder 500gr in PE bottle in colors: transparent, pink

Liquid 250 ml in glass bottle + instruction in a cupboard box



No. 33, Malekian Karimkhan Zand TEL: (+98-21) 88318502 e-mail: info@marlic.com

Alley, Iranshahr St.; Ave. Tehran- Iran Fax: (+9821) 88820843 web site: www.marlic.com