

CLINIC PRODUCTS



ACROPARS TRII

Temporary Crown & Bridge Material

Acropars TRII "cold cure acrylic for temporary crown & bridge" is based on poly ethyl methacrylate copolymers. This product is available in 4 A series color.



APPLICATION:

- 1- Making temporary crown & bridges.

ADVANTAGE:

- 1- Having no handfull.
- 2- Making low heat and no damages to teeth tissue.
- 3- Good wear resistance.
- 4- A series, 4 color.
- 5- Applicable with direct (inside mouth) or indirect (outside mouth) methods.
- 6- Normal setting time.
- 7- Excellent polishability.

PACKAGING:

Kit consists:

Powder 4 (50 g) PE bottles

Liquid 1 (100 ml) glass bottle

Glass dropper

Mixing cup

Instruction to use



No. 33, Malekian Alley, Iranshahr St.;
Karimkhan Zand Ave. Tehran- Iran
TEL: (+98-21) 88318502 Fax: (+9821) 88820843
e-mail: info@marlic.com web site: www.marlic.com