

ACROPARS TRII

Temporary Crown & Bridge Material

Acropars TRII “self cure acrylic for temporary crown & bridge” is based on poly ethyl methacrylate copolymers. This product is available in 4 A series color.

APPLICATION:

- 1- Making temporary crown & bridges.

ADVANTAGE:

- 1- Having no handful.
- 2- Making low heat and no damages to teeth tissue.
- 3- Good wear resistance.
- 4- A series, 4 color.
- 5- Applicable with direct (inside mouth) or indirect (outside mouth) methods.
- 6- Normal setting time.
- 7- Excellent polishability.

SPECIFICATION:

No	Specification	Unit	Obtained	Required
1	Pour time	Sec	115	90-120
2	Sticky time	Sec	240	180-240
3	Setting time in the thick layer	Sec	358	300-360
4	Setting time in the thick layer	Sec	450	390-450

PACKAGING:

Kit consists of:

- Powder 4 (50 g) PE bottles
- Liquid 1 (100 ml) PE glass bottle
- Glass dropper
- Mixing cup
- Instruction to use

ET73304

Marlic Medical Industries Co.

No. 33, Malekian Alley, Iranshahr
St.; Karimkhan Zand Ave.
Tehran- Iran

Tel: (+98-21) 88318502
Fax: (+98-21) 88820843
E- Mail: info@marlic.com

Zip Code: 1584786613
Web Site: www.marlic.com

