

ACROSOFT SL1

Heat Cure Acrylic Soft liner

Acrosoft SL is based on methacrylic copolymers and is designed for permanent soft relining of removable dental prostheses. Thanks to specific properties, it adheres to prosthesis strongly and all procedure of preparation and polishing are same as acrylic denture materials (Acropars 100).

APPLICATION:

- 1- For patient with sensitive tissue.
- 2- Atrophic relined ridges
- 3- For prosthesis dealing with deep under cuts of soft tissue

ADVANTAGE:

- 1- Complete adhesion to all new and old dentures.
- 2- Enough resilient to keep normal chewing ability.
- 3- It can be polished easily.
- 4- Acrosoft SL does not shrink.

SPECIFICATION:

No	Specification	Unit	Obtained	Required
1	Working time	min	12	10-14
2	Curing time	Sec	30	30
3	Softning factor	%	1.2	$0.8 \leq P < 2.5$

PACKAGING:

Kit containing:

Powder 30 gr PE bottle in transparent and pink colors.

Liquid 30 ml glass bottle

Powder measure

Liquid measure

Instruction to use

ET73306

Marlic Medical Industries Co.

No. 33, Malekian Alley, Iranshahr
St.; Karimkhan Zand Ave.
Tehran- Iran

Tel: (+98-21) 88318502
Fax: (+98-21) 88820843
E- Mail: info@marlic.com

Zip Code: 1584786613
Web Site: www.marlic.com

