

ACROSOFT TC1

Tissue Conditioner

Acrosoft TC is a tissue conditioner and functional material that restores massages and rehabilitates damaged tissues before a new denture rebase or reline. Acrosoft TC keeps soft for 15 days.

APPLICATION:

- 1- Temporary soft reline as tissue conditioner.
- 2- Functional molding material.

ADVANTAGE:

- 1- Powder and liquid will mix easily and fast,
- 2- Resulted past has a good consistency.
- 3- It does not cause any irritation.

SPECIFICATION:

| No | Specification | Unit | Obtained | Required |
|----|---------------|------|----------|----------|
| 1 | Consistency | mm | 66 | 25-75 |
| 2 | Penetration | mm | 4 | <5 |
| 3 | Sticky time | min | 4 | 4 ± 1 |

PACKAGING:

Kit consists of:

Powder 30 gr in PE bottle, transparent and pink

Liquid 30 ml in glass bottle

Powder measure

Liquid measure

Instruction to use

ET73305

Marlic Medical Industries Co.

No. 33, Malekian Alley, Iranshahr

St.; Karimkhan Zand Ave.

Tehran- Iran

Tel: (+98-21) 88318502

Fax:(+98-21) 88820843

E- Mail: info@marlic.com

Zip Code: 1584786613

Web Site: www.marlic.com

