CLINIC PRODUCTS



# GIBAGRON GZAGG

## Temporary Crown & Bridge Material

CIBACRON "cold Cure Acrylic For Temporary Crown & Bridge" is based on poly ethyl methacrylate copolymers. This product is available in 2 & 6 colors A & B series.



### **APPLICATION:**

1- Making temporary crown & bridges.



#### **ADVANTAGE:**

- 1- No harmful material gainst teeth pulp.
- 2- Low heat generation and no damages to teeth pulp.
- 3- Good wear resistance.
- 4- A and B colors (6 colors).
- 5- Applicable with direct (inside mouth) or indirect (outside mouth) methods.
- 6- Short setting time.
- 7- Excellent polishability.

#### **PACKAGING:**

Kit consists: CIBACRON C2 CIBACRON C6 Powder 2 (80 gr) PE bottles 6 (35 gr) PE

Liquid 1 (100 ml) glass bottle 1 (100 ml) glass bottle

Glass dropper Mixing cup

Instruction to use



No. 33, Malekian Alley, Iranshahr St.; Karimkhan Zand Ave. Tehran- Iran TEL: (+98-21) 88318502 Fax: (+9821) 88820843 e-mail: info@marlic.com web site: www.marlic.com