

CIBACRON C2, C6

Temporary Crown & Bridge Material

CIBACRON “Self Cure Acrylic For Temporary Crown & Bridge” is based on poly ethyl methacrylate copolymers. This product is available in 2 & 6 colors A & B series.

APPLICATION:

1- Making temporary crown & bridges.

ADVANTAGE:

- 1- No harmful material against teeth pulp.
- 2- Low heat generation and no damages to teeth pulp.
- 3- Good wear resistance.
- 4- A and B colors (6 colors).
- 5- Applicable with direct (inside mouth) or indirect (outside mouth) methods.
- 6- Short setting time.
- 7- Excellent polishability.

SPECIFICATION:

No	Specification	Unit	Obtained	Required
1	Pour time	Sec	115	90-120
2	Sticky time	Sec	240	180-240
3	Setting time in the thick layer	Sec	300	300-360
4	Setting time in the thick layer	Sec	350	390-450

PACKAGING:

CIBACRON C2 KIT:

Powder 2 (80 gr) PE bottles
 Liquid 1 (100 ml) PE glass bottle
 Glass dropper
 Mixing cup
 Instruction to use

CIBACRON C6
 6 (35 gr) PE bottles

ET73307

Marlic Medical Industries Co.

No. 33, Malekian Alley, Iranshahr
 St.; Karimkhan Zand Ave.
 Tehran- Iran

Tel: (+98-21) 88318502
 Fax: (+98-21) 88820843
 E- Mail: info@marlic.com

Zip Code: 1584786613
 Web Site: www.marlic.com

