

LABORATORY PRODUCTS



ACROPARS 100

Heat Cure Denture Acrylic

This product is designed based on methyl methacrylate copolymers, comply with ISO 1567 and ADA No. 12. ACROPARS 100 has excellent performance for preparing full and partial dental prosthesis.



APPLICATION:

- 1- Full denture prosthesis.
- 2- Partial denture.

ADVANTAGE:

- 1- Adequate time for mixing and molding.
- 2- Curing time, 30 min. in boiling water.
- 3- Low water absorption.
- 4- High mechanical properties.
- 5- Good adhesion to synthetic teeth.
- 6- Good polishing property.
- 7- Excellent color stability.
- 8- No bubble and whitening in prosthesis.

PACKAGING:

Powder 500 gr in PE bottles
 10 kg in laminate bag
Liquid 250 ml in glass bottles



No. 33, Malekian Alley, Iranshahr St.;
Karimkhan Zand Ave. Tehran- Iran
TEL: (+98-21) 88318502 Fax: (+9821) 88820843
e-mail: info@marlic.com web site: www.marlic.com