

ACROPARS Re

Self Cure Acrylic For Repair

Acropars Re is designed for repairing or extending the full or partial prosthesis and is based on methyl methacrylate copolymer. It is comply with ISO 1567 and ADA No. 12.

APPLICATION:

- 1- Repairing denture prosthesis.
- 2- Extending denture prosthesis.
- 3- Relining of denture prosthesis.

ADVANTAGE:

- 1- Enough time to repair or reline.
- 2- Very good adhetion to old prosthesis.
- 3- Good polish ability.
- 4- Good color stability.

SPECIFICATION:

No	Specification	Unit	Obtained	Required
1	Sticky time	min	5	5±2
2	Setting time	min	11.5	9±3
3	Transverse deflection	mm	3.2	2-5.5

PACKAGING:

Powder 500gr in PE bottle in colors: transparent, pink and veined pink colors Liquid 250 ml in glass bottle + instruction in a cupboard box

Tel: (+98-21) 88318502

Fax:(+98-21) 88820843

E- Mail: info@marlic.com







Zip Code: 1584786613

Web Site: www.marlic.com