

LABORATORY PRODUCTS



# ACROPARS 200

## Cold Cure Acrylic For Special Tray

This product is designed with a unique formulation to prepare special individual tray and is available in 3 colors. The tray has excellent mechanical properties and very low shrinkage.

### APPLICATION:

- 1- Special tray
- 2- Rebase



### ADVANTAGE:

- 1- Very short sticky and enough working time.
- 2- Minimal shrinkage.
- 3- Good mechanical properties.
- 4- Good impact resistance.
- 5- In different colors.

### PACKAGING:

Powder 500 gr in PE bottle; yellow, green & blue colors and  
Liquid 250 ml in glass bottle + instruction in a cupboard box  
Kit consists of:

Powder: 1000 gr in PE bottle  
Liquid: 2 glass bottles each 250 ml  
Mixing cup  
Instruction to use



No. 33, Malekian Alley, Iranshahr St.;  
Karimkhan Zand Ave. Tehran- Iran  
TEL: (+98-21) 88318502 Fax: (+9821) 88820843  
e-mail: info@marlic.com web site: www.marlic.com