CLINIC PRODUCTS



AGROSOFT TG

Tissue Conditioner

Acrosoft TC is a tissue conditioner and functional material that restores massages and rehabilitates damaged tissues before a new denture rebase or reline. Acrosoft TC keeps soft for 15 days.



APPLICATION:

- 1- Temperorary soft reline as tissue conditioner.
- 2- Functional molding material.

ADVANTAGE:

- 1- Powder and liquid will mix easily and fast,
- 2- Resulted past has a good consistency.
- 3- It does not cause any irritation.

PACKAGING:

Kit consists of:

Powder 30 gr in PE bottle, transparent and pink

Liquid 30 ml in glass bottle

Powder measure Liquid measure Instruction to use



No. 33, Malekian Karimkhan Zand TEL: (+98-21) 88318502 e-mail: info@marlic.com

Alley, Iranshahr St.; Ave. Tehran- Iran Fax: (+9821) 88820843 web site: www.marlic.com